

Core Panel Multi-Tumor Slides MULTI-TUMOR POSITIVE CONTROL TISSUE SLIDES

Technical Data Sheet

Product Category

Core Panel Multi-Tumor Positive Control Slides

Formats Supplied

- 5 slides per pack
- 25 slides per pack

PRODUCT NUMBER: TS900 (pack of 5 slides)

PRODUCT NUMBER: TS900-25 (pack of 25 slides)

~~~~~

#### **PRODUCT DESCRIPTION**

Each Multi-tumor Core panel control slide includes the following tumor tissue sections: Brain, Thyroid, Lung, Muscle, Stomach, Small Bowel, Large Bowel, Prostate, Kidney, Uterus, Liver, Pancreas, Breast and Skin.

**FIXATION and PROCESSING:** Formalin fixed, paraffin embedded.

**THICKNESS:** 4 microns

#### **APPLICATIONS and METHOD OF USE**

**Core panel multi tumor control slide** is generally employed as the **positive control and negative control on the same slide**. These multi-tumor control slides are employed in a variety of testing such as Immunohistochemistry (IHC), fluorescent labeling (IF) and in situ-hybridization (ISH). The case tissue or test sections may be also mounted on the same slide for ease of handling and microscopic examination.

**INNOVEX Control Slides as provided are NOT pre-baked in order to preserve antigenicity and to further allow the end user the option of placing case or patient tissue on these control slide if desired. The core control slide bearing case or patient tissue section can then be baked at 60°C for 1-hour.**

**Bake slides in an oven at 60°C for 1-hour and cool to room temperature for 15-minutes prior to deparaffinization. These slides can also be baked in batches weeks and months in advance.**

**When desired, the case or patient tissue section can be placed onto a different location on this core multi tissue slide and then baked.**

#### **TESTED APPLICATIONS**

IHC (paraffin), *in-situ* hybridization (ISH), Immunofluorescence (IF).

#### **STORAGE CONDITIONS**

Store at room temperature in a cool dry place.

#### **IMPORTANT NOTE:**

- **CAP Acceptable**
- The interpretation of test results is the sole responsibility of the end user.

**FOR PROFESSIONAL AND LABORATORY USE ONLY**