

# TURBO ACTION STABLE Liquid DAB and AEC for PEROXIDASE (HRP) STAINING

## ADVANCED FEATURES

- Good to the last drop, mixed left over AEC and DAB solution stable for up to 3 weeks
- Fast developing (1-5 minutes)
- Background free
- No fading with time
- Highly stable
- 2 components only
- Environmentally friendly, mixed substrate/chromogen solution reusable for up to 3 weeks
- No precipitants
- Vibrant and crisp staining, makes visualization of cell by cell possible
- Guaranteed 2 years of shelf life from delivery
- Economical



DAB staining of Pan cytokeratin antibody in 40 minutes  
with Innovex STAT-Q staining system, no digestion  
or heat treatment



SUBSTRATE/  
CHROMOGENS

## PRODUCT ORDERING INFORMATION

Product	Product #	Size	# of Slides	Shelf life
Liquid DAB/ Substrate system:	NB314SBD	70 ml	1200	3 yrs
Liquid AEC/ Substrate system :	NB314SCS	70 ml	1200	3 yrs

## FOR IMMUNO-ALKALINE PHOSPHATASE STAINING

### ADVANCED FEATURES

- Fast developing (5-10 minutes)
- Background free
- No fading with time
- Highly stable
- Duo component Innovex Fast Red
- Single solution Permanent Innovex Brown (DAB like color for alkaline phosphatase staining)
- Single solution Permanent BCIP/NBT with no background or leeching out
- No precipitants
- Vibrant and sharp staining, makes visualization of cell by cell possible.
- Guaranteed 2 years of shelf life from delivery
- Economical



Innovex Fast Red staining of PSA antibody in 50 minutes  
with Innovex HISTO-STAT Polymer alk-phos kit

### PRODUCT ORDERING INFORMATION

Product	Product #	Size	# of Slide	Shelf life
Innovex FAST RED	NB317SBF	70 ml	1000	3 yrs
Innovex Brown (single solution)	NB317B	30 ml	500	3 yrs
Permanent, Background free BCIP/NBT (single solution)	NB321NBT	10 ml	180	3 yrs

INNOVEX



SUBSTRATE/  
CHROMOGENS